REF 8319550

en	Panel for the Qualitative Detection of Drugs of Abuse in Oral Fluid
de	Gerät zum qualitativen Nachweis von Drogen im Speichel
fr	Panneau de Détection qualitative des drogues illicites dans la salive
es	Panel para la detección cualitativa de drogas en la saliva
it	Tavoletta per la determinazione qualitativa di droghe d'abuso (DOA) nella saliva
nl	Panel voor de kwalitatieve detectie van drugs in speeksel
fi	Liiska huumausaineiden kvalitatiiviseen toteamiseen syljestä
bg	Комбиниран тест за установяване употребата на наркотични вещества в проба от слюнка

Instructions for Use

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Käyttöohje

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EN – Instructions for Use

For Your Safety

Strictly Follow the Instruction for Use

Any use of the Dräger DrugCheck™ requires full understanding and strict observation of the instructions.
The Dräger DrugCheck™ is only to be used for purposes specified herein.

Liability for Proper Function or Damage

The liability for the proper function of the Dräger DrugCheck™ is irrevocably transferred to the owner or operator to the extent that the kits are used in a manner not conforming to its intended use.
Dräger Safety cannot be held responsible for damage caused by non-compliance with the recommendations given herein.
The warranty and liability provisions of the terms of sale and delivery of Dräger Safety are likewise not modified by the recommendations given above.

Dräger Safety AG & Co. KGaA

Intended Use

For Cocaine, Opiate, Amphetamine, Phencyclidine, Methamphetamine, and THC

Dräger DrugCheck™ detects drugs of abuse in oral fluid and provides for the collection, screening, transport, and storage of saliva specimens.
Dräger DrugCheck™ is intended for in-vitro diagnostic and forensic use in simultaneous, qualitative detection of multiple drugs and metabolites in oral fluid. The Dräger DrugCheck™ profile is based, when possible, on the preliminary oral fluid testing requirements of the U.S. Federal Workplace Drug Testing Programs¹⁾.

Dräger DrugCheck™ Cutoff Concentrations

Drug	ng/mL	
COC	Cocaine	20
OPI	Opiate	40
AMP	d-Amphetamine	50
MET	d-Methamphetamine	50
PCP	Phencyclidine	10
THC	Δ ⁹ -Tetrahydrocannabinol	50

Dräger DrugCheck™ provides only preliminary analytical test results. An alternate, more specific chemical method must be used to obtain confirmed analytical results. Gas Chromatography/Mass Spectrometry (GC/MS or LC/MS/MS) is the preferred confirmation method for oral fluid specimens¹⁾. Clinical consideration and professional judgment should be applied to drugs of abuse test results, particularly preliminary positive results.

Summary

Drugs may be ingested, inhaled, or injected²⁾. After entering the bloodstream, they are rapidly metabolized by various pathways³⁾. Many drugs and drug metabolites are excreted in urine and detected by traditional drug urine screening assays. Alternative laboratory methods have been used to detect drugs in blood or serum.

The utility of oral fluid as a biological specimen for the detection of recent drug use has been supported by a number of studies⁴⁾⁻⁸⁾. A growing body of literature describes the detection and related pharmacokinetics of cocaine, opiates, and amphetamines in oral fluid^{8),9)}. Correlation of drug in oral fluid and blood has been reported¹⁰⁾⁻¹¹⁾. Unlike urine drug tests, oral fluid drug testing detects active drugs present at low concentrations (recent drug use, typically within hours).

Historically, drug-screening assays relied on classical chemical methods, such as thin-layer chromatography or liquid chromatography, which are accurate but laborious procedures. More recent assays for initial testing are based on immunoassay principles. Examples of these methods are radioimmunoassay, enzyme immunoassay, kinetic interaction of microparticles in solution, and fluorescence immunoassay. These methods require sophisticated instrumentation and the handling of laboratory reagents and urine specimens. The classic drug screening technology incorporated in the Dräger DrugCheck™ provides the accuracy of immunoassays but does not require the collector to witness urine collection or handle laboratory reagents and urine specimens.

Dräger DrugCheck™ is a simple, qualitative, visually read 10-minute test that detects recent drug use using oral fluid specimens. The design permits convenient collection and lateral flow screening, with specimen available for confirmation testing.

Test Principle

Dräger DrugCheck™ assays are based on the principle of competitive inhibition. Drug that may be present in the oral fluid competes with drug conjugate immobilized on the test membrane for binding sites on the antibody-coated microparticles. A visual signal is generated through classical lateral flow technology.

Specimen Collection: Oral fluid is collected by absorption directly into the Dräger DrugCheck™ foam-tipped oral fluid collector and is delivered into the collection tube by pushing the collector, foam-first, down into the expresser. The screening test is initiated after placing the profile card into the collection tube.

The Test:

The oral fluid interacts with antibody-coated microparticles and drug conjugate present on the membrane. In the absence of drug, the antibody is free to react with the drug conjugate, causing the formation of a red/pink band in each drug zone. When drug is present in the specimen, it binds to the antibody-coated microparticles. If sufficient drug is present, the microparticles are inhibited from binding the drug conjugate, and no band is formed. A preliminary positive ("non-negative") specimen causes the detection area on the membrane to remain blank.

Quality Control:

An additional antibody/antigen reaction is designed into the Test Valid area. When the specimen application and procedure are appropriate, antibodies on the reagent membrane bind the antigen on the microparticles, and the Test Valid band forms.

Chain of Custody:

Excess oral fluid sample can be used for confirmation testing. The tube is ready for delivery to a confirmation site for laboratory testing.

Precautions & Warnings

Dräger DrugCheck™ tests are for in-vitro diagnostic use in testing human oral fluid for designated drugs.

Follow universal precautions when handling specimens and use safe laboratory procedures such as those outlined in Biosafety in Microbiological and Biomedical Laboratories (HHS Publication Number [CDC] 93-8395).

Storage & Stability

Store Dräger DrugCheck™ tests at room temperature, 15 – 30°C (59 – 86°F), in the original foil pouch. If the foil pouch containing the Dräger DrugCheck™ profile card is damaged (e.g., a hole or tear), do not use that profile card. Do not use Dräger DrugCheck™ profile card after the expiration date on the foil pouch.

Materials Provided

- 25 individually packaged Dräger DrugCheck™ Profile Cards
- 25 Oral Fluid Collectors sealed in a plastic wrap
- 25 Capped Collection Tubes with Oral Fluid Expresser

Ancillary Materials

Latex Gloves, Collectors/Expressers:

Contact your local distributor for IVD-equipment and call Dräger Safety representative or technical support for assistance.

Collect Oral Fluid & Start the Test

Collection Procedure

Prior to collection, instruct the donor to drink four (4) ounces of water, and to place nothing else in the mouth such as food, drink, gum, or tobacco products for 10 minutes. Then instruct the donor to draw a pool of oral fluid into the mouth with a "puckering" action for a few moments before collection.

Collect fresh oral fluid specimens as follows; no preservative, special handling, or pretreatment is required, and no diluents, dilution, or buffer is necessary.

1. Instruct the donor to remove the collector from the plastic wrap, and place the collector foam inside the mouth for up to three (3) minutes, completely saturating the collector foam with oral fluid.
Specifically, instruct that the donor:
 - Bathe the collector foam in the pool of oral fluid
 - Alternatively, stimulate salivation by continuously moving the collection foam along the sides and lower margin of the tongue
 - Not chew, compress, or suck on the collector foam
2. Observe the collection.
3. Have the donor remove the collector from the mouth and hand it to the administrator. The collector foam should be very wet, and must not be compressed when removed.
 - Guide the collector foam down into the oral fluid expresser at the opening of the tube, and slowly push the collector downward until it comes to a firm stop.

i **NOTE**
There will be some resistance and a popping sound about 2/3 the way down the expresser, make sure the collector is pushed past that point, all the way to the bottom of the expresser (see illustrated instructions). To express the maximum amount of sample, lift up slightly and push downward a second time.

- Oral fluid flows directly into the collection tube.
- Pull straight up on the collector and lift out of the tube; the expresser will be attached to it. Throw both of these pieces away.

i **NOTE**
It is very important that the collector remain in the donor's mouth for a sufficient time period. The goal is to collect a minimum of one (1) milliliter of sample (use the markings on the tube as a guide). If insufficient sample is collected, repeat the collection with a new collector/expresser, adding the additional sample to the same tube.

4. Open the foil pouch and remove profile card. Record donor ID if desired in the writable area and drop card down into the tube, "arrow" first. Allow testing card to drop to the bottom of the tube.

i **NOTE**
Do not place profile card in mouth.

5. Wait 10 minutes and interpret results. Do not read results after 15 minutes. See next section for interpretation guidelines.

6. If result confirmation is necessary, remove the testing card and firmly seal the tube. Place evidence tape up and over the cap, and place tube into a chain-of-custody bag. Always ship in accordance with applicable local, state, or federal regulations.

Interpret Test Validity and Results (A)

Test validity and results must be interpreted between 10 and 15 minutes after initiating the test.

1. Interpret test validity: A valid test is indicated by a band in the area of the card next to the sign "TV". The reagent strip shown on **A1** is valid. Since a valid test may give a faint or incomplete band, any Test Valid band confirms that the test is valid. Due to the high viscosity and variability of some oral fluid samples, test results may require up to 15 minutes to form.
An invalid test is indicated by the absence of a distinct band in the Test Valid area or by a reddish reagent background on the strip which obscures the presence of bands 15 minutes after the test is initiated. If an invalid result is obtained, see Step 3.
2. Immediately interpret test results as either negative or preliminary positive.
 - a. A negative result for a given drug (i.e., drug absent or below the cutoff) is the presence of a band in the area adjacent to the drug name (**A2**).

i **NOTE**
The intensity of the bands in the Results area may vary. A negative sample may give a faint or incomplete band; any band in the Result area indicates a negative result. The test results pictured on **A2** are negative for all drugs in this device.

- b. A preliminary positive result for a given drug (i.e., drug present above the cutoff, suggesting current or recent drug use) is the **absence** of a band in the test result area adjacent to the drug name. The adjacent area appears off-white. The result shown on **A3** is preliminarily positive for Cocaine.
 - c. If all drug results are negative, or if no further confirmation tests are required, properly discard the test card and tube.
 - d. When confirmation of a preliminary positive screen result is required, follow established Chain-of-Custody procedures for shipment to a laboratory.
3. If test did not run, or if invalid test results are obtained, use a second Dräger DrugCheck™ device, after discarding the first.
 - a. Instruct the donor to drink 4 ounces of water. Assure that this water is swallowed and that nothing else is placed in the mouth for three to five minutes.
 - b. During the minute or two prior to collection, instruct the donor to draw a pool of oral fluid into the mouth with a "puckering" action.
 - c. Provide a second packaged collector to the donor, and repeat collection instructions with emphasis on:
 - Bathing the collection foam in oral fluid, and/or
 - Continuous gentle movement of the collector over the tongue throughout the collection
 - d. Observe a **4-minute** collection and perform the test according to the instructions in "Collect Oral Fluid & Start Test" section.

Quality Control

Quality control testing at regular intervals is good analytical practice and may be required by federal, state, or local guidelines. Always check with the appropriate licensing or accrediting bodies to ensure your quality control program meets the established standards.

Internal control: Each Dräger DrugCheck™ profile card has a built-in process control. The Test Valid band should always appear if the test procedure is performed properly and the specimen is adequate. In addition, the red background of reagents in the read windows normally clears to off-white, providing a distinct result 10 minutes after initiating the test. A Test Valid band is an internal process control. Partial, incomplete, or broken bands still indicate a valid test.

External control: The use of positive and negative oral fluid controls is recommended to test each shipment of product, when a new lot is used, or anytime the product has been stored outside the recommended storage conditions.

Contact Dräger Safety, for information regarding commercially available controls. These controls should produce the expected results. If these results are not obtained, call Technical Support for assistance.

Limitations

1. Rare false positive or false negative results may occur as a result of nonspecific interactions (physiological variation, state of health, contamination of oral fluid with food or other substances, etc.).
2. Antibodies used in these reagents are very specific for targeted drugs of abuse. Nevertheless, structurally similar prescribed and over-the-counter drugs can react with the antibody reagents and cause false positive results. A more specific alternate chemical method such as GC/MS or LC/MS/MS must be used to obtain a confirmed analytical result. A test protocol for simultaneous GC/MS determination of these drugs in oral fluid is available e.g. from Varian, Inc., for validation and use by commercial laboratories; contact Dräger Safety, for further information. Clinical consideration and professional judgment should be applied to any drug abuse test result, particularly preliminary positive test results, and should be reviewed by experts / professionals.
3. Oral fluid collection may be closely observed. As a result, specimen adulteration is unlikely. Nevertheless, if adulteration or substituted sample is suspected, discard that specimen and collect a second specimen for testing.

Precision (T1)

The performance of the strip lots was studied over a minimum period of three days. Three operators read the results of devices run with spiked oral fluid at OX, 0.25X, 0.5X, 1X, and 1.5X of each drug cut-off in a blind, random order. All lots are required to minimally perform with a ≥ 95% confidence level that negative results will be attained with negative samples (OX synthetic), and a ≥ 95% confidence level that positive results will be attained with drugs at 150% of their respective cutoff concentrations.

Test strip lots performance results at OX and 1.5X of each drug assay cut-off are summarized on Table 1 (T1).

Specificity (T2)









Structurally similar compounds for cocaine, morphine, amphetamine, THC, methamphetamine, and phencyclidine were tested with Dräger DrugCheck™. The results in Table 2 (T2) are expressed as that amount of the compound capable of giving a result equivalent to the cutoff for that assay.

Cross Reactivity (T3)

The Compounds listed in Table 3 (T3) were tested at 100,000 ng/mL and found not to cross-react with the Dräger DrugCheck™ cocaine, morphine, amphetamine, THC, methamphetamine, and phencyclidine assays (except as noted).

References

- 1) Mandatory Guidelines for Federal Workplace Drug Testing Programs, Draft #3, 12/5/00 Meeting of the Drug Testing Advisory Board. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
- 2) O'Brien CP. Drug Addiction and Drug Abuse. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. Goodman & Gilman's The Pharmacological Basis of Therapeutics. 9th ed. McGraw-Hill Co. Inc. 1996:557-577.
- 3) Karch SB. Drug Abuse Handbook. CRC Press LLC 1998.

- Cone, E.J. Saliva Testing for Drugs of Abuse. Ann. N.Y. Acad. Sci., 694, 1993.
- Jenkins, A.J. Detecting drugs of Abuse in Saliva. Ther. Drug Monit. and Toxicol. 19(3), 1998.
- Kidwell, DA; Holland, J, et al. Testing for drugs of abuse in saliva and sweat. J. Chromatogr. B 713, 1998, 111–135.
- Peel, HW; Perrigo, BJ, et al. Detection of Drugs in Saliva of Impaired Drivers. J. Forensic Sci. 29(1), 1984, 185–189.
- Schramm,W; Smith RH, et al. Drugs of Abuse in Saliva: A Review. J. Anal. Toxicol. 16, 1992, 1–9.
- Cone, EJ; Oyler, J, et al. Cocaine Disposition in Saliva Following Intravenous, Intranasal, and Smoked Administration. J. Anal. Toxicol., Vol. 21, 1997, 465–475.
- Jenkins, AJ; Oyler, JM, et al. Comparison of Heroin and Cocaine Concentrations in Saliva with Concentrations in Blood and Plasma. J. Anal. Toxicol. 19, 1995, 359–374.
- O’Neal, CL; Crouch, DJ, et al. Correlation of Saliva Codeine Concentrations with Plasma Concentrations after Oral Codeine Administration. J. Anal. Toxicol. 23, 1999, 452–459.

Step-by-Step (1-7)

- Have donor keep foam collector in the mouth until thoroughly soaked (up to three minutes). They should not chew or suck on the foam.

- The administrator should place the collector foam-first into the expresser in the tube.

Slowly push the collector down into the expresser. You will feel some resistance about 2/3 the way down; push past this to reach the bottom of the expresser (you will hear a popping sound). The oral fluid will flow into the tube.

- To express the maximum amount of sample, lift up slightly and push downward a second time. NOTE: It is important that a minimum of 1 mL of sample is collected.

- Lift the collector out of the tube at a straight upward angle. The expresser will be removed with it. Throw both pieces away.

- Open the foil pouch and remove testing card. Record name or ID number in writable area. NOTE: Do not place profile card in mouth.

- Place testing card in tube as shown: "arrow" down. Allow testing card to drop to the bottom of the tube. Wait 10 minutes and interpret results. Do not read results after 15 minutes. The presence of any line is negative; absence of any line is positive. See "Interpret Test Validity and Results" section for guidelines.

- To send for confirmation: Remove the testing card and tightly seal the tube with the cap. Complete your facility’s Chain-of-Custody procedures and pack into a box for shipping to lab. Ship according to local, state, or federal regulations.

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Referenzen

- Mandatory Guidelines for Federal Workplace Drug Testing Programs, Draft #3, 12/5/00 Meeting of the Drug Testing Advisory Board. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
- O'Brien CP. Drug Addiction and Drug Abuse. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. Goodman & Gilman's The Pharmacological Basis of Therapeutics. 9th ed. McGraw-Hill Co. Inc. 1996:557-577.
- Karch SB. Drug Abuse Handbook. CRC Press LLC 1998.
- Cone, E.J. Saliva Testing for Drugs of Abuse. Ann. N.Y. Acad. Sci., 694, 1993.
- Jenkins, A.J. Detecting drugs of Abuse in Saliva. Ther. Drug Monit. and Toxicol. 19(3), 1998.
- Kidwell, DA; Holland, J, et al. Testing for drugs of abuse in saliva and sweat. J. Chromatogr. B 713, 1998, 111–135.
- Peel, HW; Perrigo, BJ, et al. Detection of Drugs in Saliva of Impaired Drivers. J. Forensic Sci. 29(1), 1984, 185–189.
- Schramm,W; Smith RH, et al. Drugs of Abuse in Saliva: A Review. J. Anal. Toxicol. 16, 1992, 1–9.
- Cone, EJ; Oyler, J, et al. Cocaine Disposition in Saliva Following Intravenous, Intranasal, and Smoked Administration. J. Anal. Toxicol., Vol. 21, 1997, 465–475.
- Jenkins, AJ; Oyler, JM, et al. Comparison of Heroin and Cocaine Con-centrations in Saliva with Concentrations in Blood and Plasma. J. Anal. Toxicol. 19, 1995, 359–374.
- O’Neal, CL; Crouch, DJ, et al. Correlation of Saliva Codeine Concen-trations with Plasma Concentrations after Oral Codeine Administra-tion. J. Anal. Toxicol. 23, 1999, 452–459.

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Références

- Mandatory Guidelines for Federal Workplace Drug Testing Programs, Draft #3, 12/5/00 Meeting of the Drug Testing Advisory Board. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
- O'Brien CP. Drug Addiction and Drug Abuse. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. Goodman & Gilman's The Pharmacological Basis of Therapeutics. 9th ed. McGraw-Hill Co. Inc. 1996:557-577.
- Karch SB. Drug Abuse Handbook. CRC Press LLC 1998.
- Cone, EJ. Saliva Testing for Drugs of Abuse. Ann. N.Y. Acad. Sci., 694, 1993.
- Jenkins, AJ. Detecting drugs of Abuse in Saliva. Ther. Drug Monit. and Toxicol. 19(3), 1998.
- Kidwell, DA; Holland, J, et al. Testing for drugs of abuse in saliva and sweat. J. Chromatogr. B 713, 1998, 111–135.
- Peel, HW; Perrigo, BJ, et al. Detection of Drugs in Saliva of Impaired Drivers. J. Forensic Sci. 29(1), 1984, 185–189.
- Schramm,W, Smith RH, et al. Drugs of Abuse in Saliva: A Review. J. Anal. Toxicol. 16, 1992, 1–9.
- Cone, EJ; Oyler, J, et al. Cocaine Disposition in Saliva Following Intravenous, Intranasal, and Smoked Administration. J. Anal. Toxicol., Vol. 21, 1997, 465–475.
- Jenkins, AJ; Oyler, JM, et al. Comparison of Heroin and Cocaine Concentrations in Saliva with Concentrations in Blood and Plasma. J. Anal. Toxicol. 19, 1995, 359–374.
- O'Neal, CL; Crouch, DJ, et al. Correlation of Saliva Codeine Concentrations with Plasma Concentrations after Oral Codeine Administration. J. Anal. Toxicol. 23, 1999, 452–459.

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- Mandatory Guidelines for Federal Workplace Drug Testing Programs, Draft #3, 12/5/00 Meeting of the Drug Testing Advisory Board. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
- O'Brien CP. Drug Addition and Drug Abuse. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. Goodman & Gilman's The Pharmacological Basis of Therapeutics. 9th ed. McGraw-Hill Co. Inc. 1996:557-577.
- Karch SB. Drug Abuse Handbook. CRC Press LLC 1998.
- Cone, EJ. Saliva Testing for Drugs of Abuse. Ann. N.Y. Acad. Sci., 694, 1993.
- Jenkins, AJ. Detecting drugs of Abuse in Saliva. Ther. Drug Monit. and Toxicol. 19(3), 1998.
- Kidwell, DA; Holland, J, et al. Testing for drugs of abuse in saliva and sweat. J. Chromatogr. B 713, 1998, 111–135.
- Peel, HW; Ferrigo, BJ, et al. Detection of Drugs in Saliva of Impaired Drivers. J. Forensic Sci. 29(1), 1984, 185–189.
- Schramm,W; Smith RH, et al. Drugs of Abuse in Saliva: A Review. J. Anal. Toxicol. 16, 1992, 1–9.
- Cone, EJ; Oyler, J, et al. Cocaine Disposition in Saliva Following Intravenous, Intranasal, and Smoked Administration. J. Anal. Toxicol., Vol. 21, 1997, 465–475.
- Jenkins, AJ; Oyler, JM, et al. Comparison of Heroin and Cocaine Concentrations in Saliva with Concentrations in Blood and Plasma. J. Anal. Toxicol. 19, 1995, 359–374.
- O'Neal, CL; Crouch, DJ, et al. Correlation of Saliva Codeine Concentrations with Plasma Concentrations after Oral Codeine Administration. J. Anal. Toxicol. 23, 1999, 452–459.

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Referenties

- Mandatory Guidelines for Federal Workplace Drug Testing Programs, Draft #3, 12/5/00 Meeting of the Drug Testing Advisory Board. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
- O'Brien CP. Drug Addiction and Drug Abuse. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. Goodman & Gilman's The Pharmacological Basis of Therapeutics. 9th ed. McGraw-Hill Co. Inc. 1996:557-577.
- Karch SB. Drug Abuse Handbook. CRC Press LLC 1998.
- Cone, EJ. Saliva Testing for Drugs of Abuse. Ann. N.Y. Acad. Sci., 694, 1993.
- Jenkins, AJ. Detecting drugs of Abuse in Saliva. Ther. Drug Monit. and Toxicol. 19(3), 1998.
- Kidwell, DA; Holland, J, et al. Testing for drugs of abuse in saliva and sweat. J. Chromatogr. B 713, 1998, 111– 135.
- Peel, HW; Perrigo, BJ, et al. Detection of Drugs in Saliva of Impaired Drivers. J. Forensic Sci. 29(1), 1984, 185–189.
- Schramm,W; Smith RH, et al. Drugs of Abuse in Saliva: A Review. J. Anal. Toxicol. 16, 1992, 1–9.
- Cone, EJ; Oyler, J, et al. Cocaine Disposition in Saliva Following Intra-venous, Intranasal, and Smoked Administration. J. Anal. Toxicol., Vol. 21, 1997, 465–475.
- Jenkins, AJ; Oyler, JM, et al. Comparison of Heroin and Cocaine Con- centrations in Saliva with Concentrations in Blood and Plasma. J. Anal. Toxicol. 19, 1995, 359–374.
- O'Neal, CL; Crouch, DJ, et. al. Correlation of Saliva Codeine Concen- trations with Plasma Concentrations after Oral Codeine Administra- tion. J. Anal. Toxicol. 23, 1999, 452–459.

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- Mandatory Guidelines for Federal Workplace Drug Testing Programs, Draft #3, 12/5/00 Meeting of the Drug Testing Advisory Board. Department of Health and Human Services, Substance Abuse and Mental Health Services Administration.
- O'Brien CP. Drug Addiction and Drug Abuse. In: Hardman JG, Limbird LE, Molinoff PB, Ruddon RW, Goodman LS, Gilman A, eds. Goodman & Gilman’s The Pharmacological Basis of Therapeutics. 9th ed. McGraw-Hill Co. Inc. 1996:557-577.
- Karch SB. Drug Abuse Handbook. CRC Press LLC 1998.
- Cone, EJ. Saliva Testing for Drugs of Abuse. Ann. N.Y. Acad. Sci., 694, 1993.
- Jenkins, AJ. Detecting drugs of Abuse in Saliva. Ther. Drug Monit. and Toxicol. 19(3), 1998.
- Kidwell, DA; Holland, J, et al. Testing for drugs of abuse in saliva and sweat. J. Chromatogr. B 713, 1998, 111–135.
- Peel, HW; Perrigo, BJ, et al. Detection of Drugs in Saliva of Impaired Drivers. J. Forensic Sci. 29(1), 1984, 185–189.
- Schramm,W; Smith RH, et al. Drugs of Abuse in Saliva: A Review. J. Anal. Toxicol. 16, 1992, 1–9.
- Cone, EJ; Oyler, J, et al. Cocaine Disposition in Saliva Following Intravenous, Intranasal, and Smoked Administration. J. Anal. Toxicol., Vol. 21, 1997, 465–475.
- Jenkins, AJ; Oyler, JM, et al. Comparison of Heroin and Cocaine Concentrations in Saliva with Concentrations in Blood and Plasma. J. Anal. Toxicol. 19, 1995, 359–374.
- O’Neal, CL; Crouch, DJ, et. al. Correlation of Saliva Codeine Concentrations with Plasma Concentrations after Oral Codeine Administration. J. Anal. Toxicol. 23, 1999, 452–459.

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Т1 – Precision:

d-Amphetamine	Precision (%)			
(Cut-off = 50 ng/ml)	Lot 1	Lot 2	Lot 1	Lot 2
	+	–	+	–
0	0	100	0	100
12.5	0	100	0	100
25	0	100	6.7	93.3
50	100	0	100	0
75	100	0	100	0

Cocaine	Precision (%)			
(Cut-off = 20 ng/ml)	Lot 1	Lot 2	Lot 1	Lot 2
	+	–	+	–
0	0	100	0	100
5	6.7	93.3	90	10
10	83.3	16.7	100	0
20	100	0	100	0
30	100	0	100	0

d-Methamphetamine	Precision (%)			
(Cut-off = 50 ng/ml)	Lot 1	Lot 2	Lot 1	Lot 2
	+	–	+	–
0	0	100	0	100
12.5	0	100	0	100
25	53.3	46.7	56	44
50	100	0	98	2
75	100	0	100	0

Morphine	Precision (%)			
(Cut-off = 40 ng/ml)	Lot 1	Lot 2	Lot 1	Lot 2
	+	–	+	–
0	0	100	0	100
10	0	100	6.7	93.3
20	3.3	96.7	66.7	33.3
40	86.7	13.3	100	0
60	96.7	3.3	100	0

Phencyclidine	Precision (%)			
(Cut-off = 10 ng/ml)	Lot 1	Lot 2	Lot 1	Lot 2
	+	–	+	–
0	0	100	0	100
2.5	0	100	24	76
5	6.7	93.3	94	6
10	96.7	3.3	100	0
15	100	0	100	0

Δ9-Tetrahydrocannabinol	Precision (%)			
(Cut-off = 50 ng/ml)	Lot 1	Lot 2	Lot 1	Lot 2
	+	–	+	–
0	0	100	0	100
12.5	0	100	0	100
25	0	100	6	94
50	80	20	68	32
75	100	0	96	4

T2 – Specificity:		(ng/mL)
Cocaine related compounds		
Benzoylcocgonine		50
Procaine		1000
Ecogonine		1000

Morphine related compounds	(ng/mL)
Ethylmorphine	10
6-monoacetylmorphine	10
Codeine	10
Dihydrocodeine	10
Diacetylmorphine	10
Morphine-3β-D-glucuronide	10
Hydromorphone	25
Thebaine	50
Hydrocodone	50
Oxycodone	1000
Hydroxyzine	50,000
Atropine	50,000
Nicotine	100,000
Ofloxacin	100,000

Amphetamine related compounds	(ng/mL)
d,l-Amphetamine	100

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© Dräger Safety AG & Co. KGaA